



FLOCARE® BOLUS POUCH CONNECTOR

EN • FLOCARE® BOLUS POUCH CONNECTOR. Enteral giving set for the delivery of bolus feeding. Intended to deliver tube feed formulas to a patient's enteral access device (feeding tube). The Flocare Bolus Pouch Connector connects the Nutrison Bolus Pouch with ENFit® feeding tubes. For enterally tube fed patients, specifically suitable for bolus feeding. For use by healthcare professionals or by lay users and patients after training.

WARNINGS: For enteral use only / Not for intravenous use / Choking hazard; loose parts / Single use, do not use for more than 24 hours / Do not reuse as this may create a risk of infection and cause the product not to function as it should / Do not use with an enteral feeding pump / No modification to the devices is allowed.

PRECAUTIONS: Store in a dry, clean place / Consult your healthcare professional before using the device for the first time / Notice: any serious incident that has occurred in relation to this device is to be reported to Nutricia and the local competent authority / Bolus Pouch Connector shall only be connected with ENFit enteral devices and Nutrison Bolus Pouch / Avoid any hand touch contamination throughout the procedure.

SET UP:

- (1). Wash and disinfect your hands. Please, follow local policies for PPE (Personal Protective Equipment)
- (2). Close the clamp.
- (3). Remove the ENPlus protective cap and dispose, securely connect the ENPlus pouch connector to the Nutrison Bolus Pouch
- (4). Remove the ENFit protective cap and dispose
- (5). Open clamp
- (6). Fill the Bolus Pouch Connector by gently squeezing the Nutrison Bolus Pouch. Close the clamp again

USE:

- (7). Connect the Bolus Pouch Connector to the enteral feeding tube
 - (8). Open the clamp on the Bolus Pouch Connector and squeeze the pouch with steady, consistent pressure to deliver nutrition. Always monitor the patient during feeding.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local policy.

AUS: Nutricia Australia Pty Ltd, The Zenith, Tower B, Level 12, 821 Pacific Highway, CHATSWOOD, NSW, 2067 AUSTRALIA 1800 060 051, www.nutriciamedical.com.au
NZL: Danone Nutricia NZ Ltd, 56-58Aintree Ave, Airport Oaks, Auckland 2022, New Zealand 0800 636 228, www.nutriciamedical.co.nz.
GBR: Nutricia Ltd, White Horse Business Park, TROWBRIDGE, Wiltshire, BA140XQ, England. Tel: 01225 711 688, www.nutricia.co.uk Patient/Carer Information Number (UK only) Tel: 08457 623 653.
IRL: Nutricia Ireland Ltd., Block 1, Deansgrange Business Park, Deansgrange, Co. DUBLIN, Ireland. www.nutricia.ie



Manufacturer: Nutricia Medical Devices BV, Taurusavenue 167, 2132 LS Hoofddorp, The Netherlands, www.nutriciaflocare.com
Production location: Nutricia Pharmaceutical Wuxi Co., Ltd., P.R. China



August 2023
200810_A