EN • FLOCARE® ENTERAL FEEDING CONTAINER. Empty enteral feeding reservoir for administration of enteral nutrition or liquids via an enteral giving set into the patient gastrointestinal tract.

WARNINGS: For enteral use only / Choking hazard; loose parts / Non-sterile / Single patient use, do not use for more than 24 hours / Do not put in microwave / Connect the container with enteral feeding devices only.

PRECAUTIONS: Store in a dry, clean place / Consult your healthcare professional before using the device / Notice: any serious incident that has occurred in relation to this device is to be reported to Nutricia and the local competent authority.

SET UP:

- (1). Remove the lilac cap.
- (2). Ensure the container is placed on a stable, level surface before filling. Fill the container carefully to avoid instability or tipping. Maximum temperature of the fluid should not exceed 60 degrees.
- (3). Volume (ml) indicated on the container is approximate.
- (4). Firmly screw on the lilac cap and remove the small cap.
- (5). Firmly screw the giving set to the container.
- (6). Turn the container upside-down and hang it. Consult your healthcare practitioner for the maximum hanging time.

USE:

Regularly check the feeding container to ensure the flow rate is correct. Always monitor the patient during feeding. After use, dispose the product and packaging in accordance with hospital, administrative and/or local policy.

GBR: Nutricia Ltd., White Horse Business Park, TROWBRIDGE, Wiltshire, BA14 0XQ,
England. Tel: 01225 711 688, www.nutricia.co.uk
Patient/Carer Information Number (UK only) Tel: 08457 623 653
IRL: Nutricia Ltd., Block 1, Deansgrange Business Park, Deansgrange, Co. Dublin, Ireland. www.nutricia.co.ie
Danone Nutricia NZ Ltd, 56-58 Aintree Avenue, Airport Oaks, Auckland 2022, New Zealand.
New Zealand Tel: 0800 636 228. www.nutricia.co.nz
Australia Tel: 1800 060 051. www.nutricia.com.au

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FIN: Nutricia, c/o Danone Oy, Siltasaarenkatu 8-10, 00530 Helsinki. Puh. (02) 274 4111. www.nutricia.fi